# 





# Difficulties related to large restorations



Complex structure of dental tissues, including enamel, dentin, cementum and pulp, is characterized by unique mechanical properties.

Clinical success of large-scale restorations depends mostly on the structure of remaining tooth and the choice of materials and their placement to mimic the original structure.

# Large-scale restorations

## What is the main reason for failure?



Load bearing capacity of remaining tissues, adhesive joint or the filling material is exceeded due to stress concentration.

## **Biomechanics of large-scale restoration failure**

Large restoration is prone to stress concentration and subsequent failure



## **Biomechanics of large-scale restoration failure**





Stress concentration zones in the cavity

concentrated in the proximal boxes" (Ausiello P. et al.: Mechanical behavior of bulk direct composite versus block composite and





## **Biomechanics of large-scale restoration failure**







NAT

2,7



## Cusp deflection

MO

MOD

SLOT

Pascal Magne: CT Scan–Based Finite Element Analysis of Premolar Cuspal Deflection Following Operative Procedures, The International Journal of Periodontics & Restorative Dentistry, Volume 29, Number 4, 2009



# The natural approach



The tooth responds to stress as a graded material, DEJ distributes the stress across a larger zone and prevents the development and propagation of cracks to dentin.

> DEJ is a connection of two different tissues: hard and brittle enamel and soft and resilient dentin.

# **Biomimetic solution**

## Reinforcing membrane distributes stress across larger zone



## Stress distribution without membrane

Stress distribution with membrane

## Fibrafil® CUBE



- The only resin composite with integrated continuous membrane which distributes stress and mimics the function of dentinoenamel junction (DEJ).
- Extends possibilities for large restorations and build-ups of endodontically treated teeth with severe loss of hard dental tissues.
- Ideal for dentin replacement in large posterior applications (pre-endodontic and postendodontic build-ups and fillings).



# Fibrafill® CUBE



## Individually-packed application doses



# Adaptation of Fibrafill® CUBE



## Easy adaptation of Fibrafill<sup>®</sup> CUBE





# Fibrafill<sup>®</sup> CUBE application principles

Integrated glass fiber membrane

Adaptation of the first Fibrafill CUBE dose



Adhesive joint

Integrated glass fiber membrane

> Adaptation of the second Fibrafill CUBE dose

Layering of Fibrafill® CUBE units with reinforcing membrane – distributing the stresses along the bottom and walls of the cavity



# Fibrafill<sup>®</sup> CUBE application principles

Incremental adaptation of individual Fibrafill® CUBE units along the bottom-wall curvature of the cavity





# Fibrafill<sup>®</sup> CUBE application principles



Magnitude and distribution of stress [MPa] in model cavity (FEA, loading by steel ball  $\emptyset$  2 mm). Integrated membrane helps to dissipate stress generated by static loading and reduces the shear stress in the adhesive joint, improving the marginal integrity of the filling. Associate Professor Zdenek Horak, Ph.D., the college of polytechnics, Jihlava, Czech Republic

Clinical protocol by MDDr. Tomas Slavicek, Brno, Czech Republic



Adaptation of Fibrafill® CUBE along the bottom-wall curvature of the cavity.

## Case report – Fibrafill® CUBE

Courtesy MDDr. Tomáš Slavíček, Brno, Czech Republic



## Restoration of large MOD cavity.

Reconstruction of approximal

wall (transformation to first-

class cavity)

Clinical protocol by MDDr. Tomas Slavicek, Brno, Czech Republic

Initial situation (old amalgam filling removal)

MDDr. TOMAS SLAVICE CMDD

## Case report – Fibrafill® CUBE

Fibrafill® CUBE units adapted along the bottom-wall curvature of

**Final restoration** 



Courtesy MDDr. Tomáš Slavíček,

Restoration of large MO cavity.



Clinical protocol by MDDr. Tomas Slavicek, Brno, Czech Republic

Initial situation (old filling removed due to the progression of secondary caries)

Reconstruction of approximal wall (transformation to firstclass cavity)





# Case report – Fibrafill® CUBE

Insertion and adaptation of Fibrafill® CUBE units along the bottom-wall curvature of the cavity.

**Final restoration** 



Courtesy MDDr. Tomáš Slavíček, Brno, Czech Republic

Restoration of large MO cavity.



Initial situation (endodontically treated tooth)

Moisture management, cleaning of cavity, application of flowable liner etc.





# Case report – Fibrafill® CUBE

## Clinical protocol by MDDr. Tomas Slavicek, Brno, Czech Republic

Application and adaptation of Fibrafill® CUBE (support of weakened walls)

Insertion of FRC post, followed by prosthetics



Brno, Czech Republic

## Post-endo reconstruction





Clinical protocol by MDDr. Tomas Slavicek, Brno, Czech Republic

Initial situation

Excavation of carious lesions, application of matrix

**Reconstruction of** approximal wall



## Case report – Fibrafill<sup>®</sup> CUBE

Application and adaptation of Fibrafill® CUBE (support of reconstructed approximal)

Insertion of FRC post, followed by prosthetics

Courtesy MDDr. Tomáš Slavíček, Brno, Czech Republic

Post-endo reconstruction





## Clinical protocol by MDDr. Tomas Slavicek, Brno, Czech Republic

## Initial situation

## Excavation of carious lesions, cleaning, application of circular matrix







## Case report – Fibrafill<sup>®</sup> CUBE

Appplication of flowable liner, insertion and adaptation of Fibrafill® CUBE along the bottom-wall curvature

Finishing by flowable composite, RCTR and indirect prosthetics

Courtesy MDDr. Tomáš Slavíček, Brno, Czech Republic

Pre-endo reconstruction







## Case report – Fibrafill® CUBE

## **RTG** contrast

## The full package includes 36 application doses (cubes) in 9 blisters



# Fibrafill® CUBE



- Distribution of stress, integrated membrane prevents stress concentration (mimics the function of dentinoenamel junction of intact tooth).
- Reduced risk of crack development and propagation through the interface between restorative material and remaining dental tissues.
- Reduced risk of severe failure of remaining hard dental tissues.
- Increased fracture toughness of the filling or build-up.
- Easy handling/packability due to optimized viscosity.
- Reduced risk of material contamination due to discrete dosing of individual units.
- Wide range of applications, economy of treatment.

## Clinical arguments:





## Cavity preparation, etching, bonding

## Fibrafill® CUBE: workflow



## Fibrafill® CUBE: workflow





## Manipulation with application doses





Adaption of individual doses in the cavity. The number of doses used depends on the size of the cavity. Application of layer of a standard composite (the composite should create 1–2 mm thick layer on the occlusal surface).

## Fibrafill® CUBE: workflow



# Fibrafill® CUBE: workflow







# Coming soon...

## Fibrafill<sup>®</sup> DENTIN

## Condensable short fibre-reinforced composite (SFRC)



## Fibrafill<sup>®</sup> FLOW

## Flowable short fibre-reinforced composite (SFRC)

